Instructions for use for: **SP844, CPJ840J6 and SP840**

Physiological Pressure Transducer with Dome and Monitoring Line Set

Intended purpose: The Physiological pressure transducer, dome and monitoring line set are intended for measurement of Invasive Blood Pressure.

-  Device catalogue number reference
-  Medical Device
-  Device Serial Number
-  Date of manufacture is found on the label
-  Indicates the need for the user to consult the instructions for use.
-  Distributor info
-  Unique device identification UDI Carrier
-  Store the Transducer at -25 - 70°C. For domes and line set see page 3.
-  Caution – placed near item that needs operator awareness

 Memscap AS
Langmyra 11
3185 Skoppum
Norway
Tel: +47 33 08 40 00

 memscap.com
e-mail: sensors@memscap.com

IFU status IFU33259 Rev00. 2025.09.30

Clinical Benefits

The clinical benefit of using our invasive blood pressure monitoring device during surgery or long-term surveillance is the accurate blood pressure monitoring that also provides the continuous wave of the pressure to be presented to the health care professionals on a monitor. During surgery or long-term surveillance, this information for instance allows the health care professionals to react to changes in a timely manner, makes the user able to maintain the patient’s blood pressure at the correct level or to confirm if the blood pressure is expected to fluctuate rapidly or significantly. Thus, our device could provide the means to save lives.

Monitoring Line Set and Domes

The physiological pressure transducer SP844, CPJ840J6 or SP840 shall be used with standard disposable dome with its patented solution against the transducer. The dome is part of the monitoring line sets also specified below. The dome is designed to fit the tubing and to be integrated into the applicable monitoring line set. Since dome in combination with monitoring line set will affect the overall performance, only accessories by model designation shall be used.

Transducer: CPJ840J6 / SP840	Transducer: SP844
Dome Model: CPJ84022, 840-22 Line Set Model: CPJ84026	Dome Model: 844-28, 844-31, 844-33, 844-34 Line Set Model: 844-26

Standards

EU Declaration of conformity contains references to harmonized standards and CS.

Single registration number

EU Declaration of conformity contains the SRN.

Intended patient population

The transducer, monitoring line set and dome are intended for adults.

Contraindications

The devices are not intended for use on other patients than adults.

The devices are not intended for transfusion of blood, blood products or medication.

Precautions and warnings for Transducer 

The intended user for these devices is qualified personnel in professional healthcare facility environments. Any serious incident to the patient in relation to these devices should be reported Memscap and the competent authority of the Member State in which the user and the patient is established.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the transducer. Otherwise, degradation of the performance of this equipment could result.

WARNING: A change of cables or connectors may affect the EMC for the transducer. Degradation of the performance of this equipment could result.

DO NOT Use cleaning agents containing iodide.

DO NOT Autoclave or use an ultrasonic process on the transducer. Maximum temperature 70°C.

DO NOT Immerse connector into any kind of liquid.

DO NOT Pull cable, either at transducer or connector end.

DO NOT Expose the metal diaphragm of the transducer to any hard, sharp or pointed objects, or to unnecessary forces. When the transducer is not being used, the protection cap should always be mounted.

MAKE SURE

- That the transducer is zeroed and calibrated on the monitor.
- **For the combination with dome or Line-set, MAKE SURE:**
- That the dome has not been pressurized before it is mounted on the transducer.
- That if a syringe is used to flush the dome and the tubing system, it should be 10 ml or larger.
- That the port to atmosphere is open on the transducer dome when flushing the dome. Flushing slowly will remove bubbles more efficiently.
- That flush fluid at room temperature or warmer is used, as warm flush fluid will decrease the occurrence of spontaneous bubbles.

- That all bubbles have been flushed completely from the monitoring line set before use.
- That the dome or monitoring line set is correctly attached to the transducer.

PLEASE NOTE

- Electrical failures may result in zero shift, incorrect signal or loss of signal.
- Use only MEMSCAP specified adapter cable, when required.

Residual risk

The transducer is optimized to ensure acceptable safety in terms of use and contamination. User attentiveness is still essential for safe operation of the transducer, as indicated by warnings and cautions in this IFU. Special attention from users is required in the following situations:

Risks related to electrical failure

Un-noticed damage or internal component failure could lead to transducer being unable to zero or failing to provide a signal.

Risks related to mechanical issues

Damage to the transducer diaphragm or incorrect connection between dome and transducer may cause incorrect signal output.

Risks related to improper use

Not adjusting the height of the monitoring line set to the same level as of the organ may cause inaccurate signal. Zero adjustment with incorrect pressure may cause incorrect signal.

Risks related to Monitor and catheter.

Please refer to the specific IFU for the monitor, saline bag, saline solution and catheter as these are not applied nor controlled by Memscap. The risks and residual risks connected to these components are stated on the respective medical device IFU.

START UP PROCEDURE

1. Turn on the monitor and connect the transducer to the pressure registration unit. Make sure that the connector is fitted correctly.
2. Activate, and follow the start-up procedure for the monitor.
3. Remove the protection cap from the transducer and slide it down the cable.
4. Take the monitoring line set out of the sterile, undamaged package and check that the configuration of the system components are correct. Ensure that all connections are safely connected and tight. Connect the monitoring line set to the transducer.
5. Close the roller clamp (2). Push the spike of the drip chamber (1) into the normal saline solution container with a twisting motion and fill the drip chamber approximately 1/3 full.
6. Place the pressure cuff around the saline container. Gravity alone is normally enough to fill the monitoring system.

7. Make sure that the roller clamp (2) and the stopcocks (6) and (9) are open. Remove the protection cap from the luer lock attachment cone on the patient tube (11).

Engage the flush rate by pressing the clip on the flush device (4). Fill the flushing system, the dome (5), the stopcocks (6) and (9) and the patient tube gently.

To minimize the formation of bubbles, use fluid at room temperature or above and fill the system slowly. When this procedure is correctly followed, the tubing system assembled and filled should be free of bubbles, and the dome diaphragm should not be distended.

8. Mount the monitoring line set with the dome on the transducer (keep the filled line set vented to atmospheric pressure).
9. Mount the transducer at the level where the pressure is to be measured.
10. Zero and calibrate the monitor as per the manufacturer's instruction.

Cleaning instructions

After use, the transducer shall be cleaned by wiping the exterior surfaces with a soft, lint-free cloth slightly dampened with clean water or a mild detergent solution (Compatible cleaning agents: Water, Isopropanol, Super Sani-Cloth Wipes). Always rinse the transducer in water after the use of any cleaning agents. The transducer shall be dried, and stored with the protective cap in place, preferably in its original packaging.

Waste management

If you need to throw away the transducer, it shall be treated as electrical and electronic equipment.

General Specifications	
Parameter	Range
Pressure Range (FSO)	-30 to 300 mmHg
Over Pressure (Max)	10 000 mmHg
Sensitivity	5 µV/V/mmHg, ±3%
Electrical excitation	4-8 VDC or -AC up to 5 KHz
Input Resistance (Min)	200 Ω
Output Resistance (Max)	3000 Ω
Non-linearity (Max)	0,5% of FSO
Zero Balance (Max)	±50 mmHg
Thermal Sensitivity Shift (Max)	0,15%/°C
Thermal Zero Shift (Max)	0,25 mmHg/ °C
Operating temperature	+15 to 40°C
Storage temperature	-25 to +70°C
Isolation Resistance (Min)	1GΩ
Leak Current (Max)	1,5 µA @250V, 50Hz
Weight	28 gram (without cable)


For CPJ840J6, CPJ84022 and CPJ84026:
 To order Philips Healthcare supplies online:
<https://www.patientcare.shop.philips.com>

Monitoring Line Set 844-26 / CPJ84026 and domes specific IFU

Additional explanation of symbols

-  Manufacturer's batch code is found on the label.
-  The monitoring line set is for single use only.
-  The expiration date is found on the label.
-  Store the package at 18-25°C.
-  Keep away from sunlight.
-  Keep dry.
-  Do not use if the package is damaged or unintentionally opened before use.
-  Single Sterile barrier system. Sterilized using Ethylene Oxide.
-  Do not re-sterilize
-  A number is inserted in the symbol to indicate the number of pieces in the package.

Position	Part description
1	Drip Chamber
2	Roller Clamp
3, 7, 8, 10	Luer Lock attachment cone
4	Flush Device
5	Dome CPJ84022 in CPJ84026 or 844-28 in 844-26
6, 9	Stop Cock
11	Luer Lock attachment Cone with Cap

General

The monitoring line set is used for invasive blood pressure monitoring and measurements in combination with the reusable physiological pressure transducer. The monitoring line set is intended for contact duration >24 hours < 30 days (Prolonged).

Cleaning instructions

The monitoring line set and dome are single-use only and shall be disposed of in accordance with the waste management procedures

Waste management

The monitoring line set and accompanying dome shall be segregated as Infectious Waste after use.

Precautions and warnings for Monitoring line set and dome

WARNING: Domes and monitoring line sets are single use only.

WARNING: Do not re-sterilize the monitoring line set or dome.

For more precautions refer to precautions for transducer where the precautions for the combination are listed under **MAKE SURE** on page 1.

Residual risk

The monitoring line set is optimized to ensure acceptable safety in terms of use and contamination. User attentiveness is still essential for safe operation of the line set, as indicated by warnings and cautions in this IFU. Special attention from users is required in the following situations:

Risks associated with damaged packaging or content

The monitoring line set could be damaged if stored incorrectly. This could lead to degradation of material causing inaccurate reading.

The monitoring line set may not remain sterile after its shelf life has expired.

The monitoring line set packaging could be damaged during packaging, transport or harsh handling. This could lead to breach of sterile barrier causing contamination, infection and disease.

Risks related to incorrect filling and flush

Flush with too little volume, rapid flush, flush with port to atmosphere closed, or flush with cold water is associated with risk of bubbles in line set.

Incorrect filling and flush may cause damage to the silicone diaphragm resulting in leakage or error reading.

Risks related to re-use of monitoring line set or dome

Re-use of the monitoring line set or dome may cause cross contamination and may result in subsequent health deterioration or even death of patient.

Risks related to improper use

Kinks in the tube/lines may cause incorrect signal.

Risks related technical issues

Un-noticed damage or internal component failure could lead to system being unable to zero, failing to provide a signal, or providing incorrect signal.

Risks related to mechanical issues

Separation of components in monitoring line set during use and may lead to risk of patient harm.

Risks related to use of other line set than specified above.

If using a different monitoring line set than CPJ84026 or 844-26, Memscap are not the owner of the risks and residual risks with this line set.

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